

Bluetooth Enabled Pharmacy Control Plus Range

400L Bluetooth Pharmacy
Freestanding Fridge with Solid Door

GDPA code PPSR400BT

SKU code 444411317



Key features

- ✓ Local Monitoring and Alarm History via Bluetooth
- ✓ Power Failure Alarm with battery back up. Rechargeable battery
- ✓ Solid and glass door options
- ✓ Intelligent Fan Management for Enhanced Performance



New Controller

Offers Local temperature monitoring and alarm history via Bluetooth. Power failure alarm with battery backup, rechargeable battery and intelligent fan management for enhanced performance



Flexible internal storage

Multiple shelf locations along with a wire basket provides a flexible storage solution.



Anti-microbial handle

Powder coated to help prevent cross contamination by inhibiting the growth of microbes.



Dual probe

Greater temperature accuracy given from two temperature probe readings. One probe measuring air temperature and one measuring load (or product) temperature.

Specification: Upright 400L Bluetooth Pharmacy Freestanding Fridge with Solid Door

Gross Capacity (litres)	394L
Net Capacity (litres)	389L
Unit External Width (mm)	595mm
Unit External Depth (mm)	660mm
Unit External Height (mm)	1935mm
Product Weight (kg)	61.8kg
Packed Width (mm)	640mm
Packed Depth (mm)	660mm
Packed Height (mm)	1995mm
Packed Weight (kg)	70.8kg
Country of Origin	UK
Internal Dimensions (mm) (WxDxH)	515mm W x 480mm D x 1783mm H
Temperature Range	2°C to 8°C
Colour	White/green
External Finish	Painted Steel
Internal Finish	HIPS*
Type of plastic	SRL800 (HIPS/ABS mix)
Door Finish	Painted Steel / Solid
Door Type	Solid
Door Style	Hinged, reversible
Door Liner	Flat
Door Open Depth (mm)	1150mm
Replacement Door Seal	Yes
Open Door Sensor	Yes
Lock	Yes, 2 Keys
Adjustable Feet (number & position)	2 adjustable (front)
Porthole	Yes
Porthole Size	16 mm
Porthole Location	Rear Bottom
Temperature Display on Panel	Yes
Temperature Alarm	Yes
Remote Alarm Contacts	Yes
Audible Alarms	Yes
Visual Alarms	Yes
Number of Shelves	7
Number of Shelf Positions	12
Shelf Size (mm) (WxDxH)	496mm W x 405mm D x 6mm H
Shelf Material	Polythene Coated Wire
Number of Baskets	1
Basket Material	Polythene Coated Wire
Basket Size (mm) (WxDxH)	485mm W x 255mm D x 171mm H
Interior Light	Yes
Light Type	LED
Light Activation	Door Swith / Controller
Fan	Yes - Tangential
Spark Free	No
Defrost setting	Auto (Off Cycle)
Controller Type	Digital BT, external
Number of Probes	2
Probe Type	Air & Load (Silicone Oil)
Ambient Temperature Range	10°C to 32°C
Refrigerant Type	R600a
CFC Free	Yes
HFC Free	Yes
HCFC Free	Yes
Ammonia Free	Yes
UK 3 Pin Plug	No
Other Plug Types Available	Yes
Voltage Supply	220 -240V
Insulation Thickness (mm)	45mm
Insulation Type	High Density Polyurethane Foam
Blowing Agent	Cyclopentane
Anti Microbial Handle	Yes
Wall Brackets	No
Optional Drawer System	No
Optional Digital Lock	No
Other Door Option	Glass PPSR394BT
Battery Backup	Yes



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.