



Eco-Drying Cabinets

ENERGY-EFFICIENT DRYING WITHOUT COMPROMISE



Tackling Sustainability with our new Drying Cabinet!

Through recent innovation, our new Eco-Drying Cabinets can significantly reduce power usage compared to standard glassware drying cabinets.

Complete with digital temperature controller and electronic 7-day timer, the eco-drying cabinet range is designed to offer market leading operational efficiencies over standard drying cabinets while remaining economically priced.

Features Include:

- Energy saving of up to 65% when compared to standard glassware drying cabinets
- Insulated on 5 sides
- Low heat output
- Digital PID temperature control
- Programmable 7-day timer

TECHNICAL INFORMATION

Cat. No	Capacity (Litres)	Air Circulation	Max Temp (°C)	Internal Dims HWD (mm)	External Dims HWD (mm)	Shelves/ Positions
TD0100SN	100	Natural Convection	85	440 x 775 x 300	680 x 900 x 390	2/20
TD0180SN	180	Natural Convection	85	555 x 995 x 325	775 x 1125 x 420	2/34

TYPICAL ENERGY USAGE AT 75 (°C)

Cat. No	Energy Used in 1hr (kWh)	Energy Used in 8hr (kWh)	Energy Used in 24 hr (kWh)
TD0100SN	0.201	1.608	4.824
TD0180SN	0.312	2.496	7.488

HEAT UP TIMES & ENERGY USAGE (TD0100SN)

Set Temp	Load	Energy Used (kWh)	Time taken/ mins
22-75°C	No	0.193	20

Tests carried out at an ambient temperature of 22°C



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.