



## RLHE0845 -40°C Sparkfree Chest Freezer

**Warning: Do not store open containers of volatile substances in this freezer.**

Compact chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

| Part Number                 | RLHE0845           |
|-----------------------------|--------------------|
| Capacity (Litres)           | 237                |
| Style                       | Chest Freezer      |
| Exterior Dimensions (HxWxD) | 860 x 1050 x 730mm |
| Interior Dimensions (HxWxD) | 640 x 845 x 445mm  |
| Weight                      | 120kg              |
| Operating Temperature       | -20°C to -45°C     |
| Interior Construction       | Aluminium          |
| Exterior Construction       | Painted Steel      |
| Defrost                     | Manual             |
| Sparkfree Interior          | ✓                  |
| High Temperature Alarm      | ✓                  |
| Low Temperature Alarm       | ✓                  |
| Digital Temperature Display | ✓                  |
| Door Lock                   | ✓                  |
| Feet                        | 4                  |
| Energy Consumption          | 0.61 kWh/24h       |
| Eco Friendly Refrigerant    | ✓                  |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.