



Warning: Do not store open containers of volatile substances in this freezer.

Compact chest freezer featuring a spark free interior for laboratory safety.

Fitted with a lid lock as standard, this freezer also comes complete with a digital temperature display and high and low temperature alarms.



Bespoke racking can be provided if required as can baskets at an extra cost. To increase the life of your freezer and to ensure you get optimum performance, preventative maintenance contracts are also available. Why not ask for details when you order?



Digital Temperature Display

RLHE1145 -40°C Sparkfree Chest Freezer

Part Number	RLHE1145
Capacity (Litres)	314
Style	Chest Freezer
Exterior Dimensions (HxWxD)	860 x 1300 x 730mm
Interior Dimensions (HxWxD)	665 x 1100 x 450mm
Weight	58kg
Operating Temperature	-20°C to -45°C
Interior Construction	Aluminium
Exterior Construction	Painted Steel
Defrost	Manual
Sparkfree Interior	✓
Door Open Alarm	✗
High Temperature Alarm	✓
Low Temperature Alarm	✓
Digital Temperature Display	✓
Key Operated Lid Lock	✓
Feet	4
Energy Consumption	0.71 kWh/24h
Eco Friendly Refrigerant	✓



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.