



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

#### Advanced Pharmacy Fridges

## Advanced Pharmacy Fridge RPFRO5044

Robustly constructed from stainless steel this refrigerator from the Labcold Advanced pharmacy range is designed to provide close temperature control for temperature sensitive vaccines, medicines and reagents.

The fan circulated interior ensures that there are no hot or cold spots and the robust construction means it's ideal for busy pharmacies where the door is frequently opened. Fitted with a sturdy door lock for extra security, remote alarm contacts and an access port are also fitted as standard.

For extra peace of mind, this Advanced Pharmacy refrigerator comes complete with a calibration certificate from our UKAS temperature calibration laboratory #8898 and is backed by a 2 year parts and labour warranty.\*

\*Parts only warranty outside of the UK

### Specifications

Product Code	RPFRO5044
Gross Capacity	150 litres
Style	Undercounter
UK Warranty	2 years
External Dimensions (mm) HxWxD	835 x 600 x 600
Temperature Range	2°C to 8°C
Fridge Set Point	5°C
Internal Material	Stainless Steel

## Specifications (cont.)


External Material	Stainless Steel
Shelves	2
Shelf Dimension (mm) WxD	475 x 385 475 x 245
Door Style	Solid
Lock Type	Key lock
Net Weight	49kg
Digital Temperature Display	Yes
High/Low Temperature Alarm	Yes
Door Open Alarm	Yes
Access Porthole	Yes
Access Porthole Size	12mm
Fan Assisted Cooling	Yes
Internal Light	Yes
Remote Alarm Contacts	Yes
Auto Defrost	Yes
Free UKAS accredited temperature calibration Certificate	Included
Power Consumption Kwh/24hrs	1.8
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

## Features

- Stainless steel construction
- High temperature alarm
- Door open alarm
- Fan assisted cooling
- Access port
- Free calibration certificate from our UKAS accredited calibration laboratory #8898
- Digital temperature display
- Low temperature alarm
- Door lock
- Internal LED light
- Remote alarm contacts
- 2 year UK parts and labour warranty

## Models

---

	Model	Gross Capacity	Door Style	Lock Type
	Advanced Pharmacy Fridge RPFRO5044	150 litres	Solid	Key lock



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.