



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

#### Advanced Pharmacy Fridges

## Advanced Pharmacy Fridge RPF21043

The white and blue enamelled exterior is robust enough to withstand frequent use and the unit comes complete with anti – tilt stainless steel shelves capable of holding up to 70kg, making these medicine fridges ideal for a busy hospital or wholesale pharmacy.

Heavy duty castors mean the Labcold RPF21043 Advanced Pharmacy Fridge is extremely manoeuvrable and easy to clean behind. In addition, the stainless steel interior is hardwearing and hygienic.

This heavy duty refrigerator can be fitted with an access port and remote alarm contacts at time of manufacture.

### Specifications

Product Code	RPF21043
Gross Capacity	650 litres
Style	Free standing
UK Warranty	2 years
External Dimensions (mm) HxWxD	1990 x 710 x 800
Temperature Range	2°C to 8°C
Fridge Set Point	4°C
Internal Material	Stainless Steel

## Specifications (cont.)





External Material	White Polyester Coated Sheet Steel
Shelves	4
Shelf Dimension (mm) WxD	500 x 620
Door Style	Glass
Lock Type	Key lock
Net Weight	138kg
Digital Temperature Display	Yes
Digital Temperature Control	Yes
High/Low Temperature Alarm	Yes
Min/Max Temperature Recording	Yes
Door Open Alarm	Yes
Access Porthole	Optional Extra
Fan Assisted Cooling	Yes
Internal Light	Yes
Lockable Castors	Yes
Remote Alarm Contacts	Optional Extra
Auto Defrost	Yes
Free UKAS accredited temperature calibration Certificate	Included
Power Consumption Kwh/24hrs	2.43
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

## Features

- Precise microprocessor temperature control
- High/low temperature alarm
- Automatic defrost
- Lockable castors
- Self closing door
- 2 year parts and labour warranty (UK only)
- Digital temperature display
- Door open alarm
- Adjustable shelving
- Door lock
- Free calibration certificate from our UKAS accredited calibration laboratory #8898

## Models

---

	Model	Gross Capacity	Door Style	Lock Type
	Advanced Pharmacy Fridge RPF21043	650 litres	Solid	Key lock
	Advanced Pharmacy Fridge RPF21043	650 litres	Glass	Key lock
	Advanced Pharmacy Fridge RPF44043	1300 litres	Solid	2 key locks
	Advanced Pharmacy Fridge RPF44043	1300 litres	Glass	2 key locks



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.