



- Extensive Range
- 90+ Years of History
- Extended Aftercare Support
- UKAS Laboratory

Advanced Pharmacy Fridges

Advanced Pharmacy Fridge RPF44043

Large double door advanced pharmacy refrigerator designed for the storage of large quantities of temperature sensitive pharmaceuticals between 2°C and 8°C.

Featuring microprocessor temperature control and ducted air flow, you can be confident that despite the size of the fridge, everything within it is stored in optimum conditions.

To assist with cleaning this Labcold vaccine fridge is fitted with heavy duty lockable casters and the interior is constructed from hygienic stainless steel.

Robustly constructed with anti tilt shelves designed to take up to 70kg, this refrigerator can withstand the rigours of any busy hospital or wholesale pharmacy.

Can be fitted with an access port and remote alarm contacts if specified at time of manufacture and come complete with a manufacturer's calibration certificate

Specifications

Product Code	RPF44043
Gross Capacity	1300 litres
Style	Free standing
UK Warranty	2 years
External Dimensions (mm) HxWxD	1990 x 1420 x 800
Temperature Range	2°C to 8°C
Fridge Set Point	5°C
Internal Material	Stainless Steel





Specifications (cont.)

External Material	White Polyester Coated Sheet Steel
Shelves	8
Shelf Dimension (mm) WxD	500 x 620
Door Style	Glass
Lock Type	2 key locks
Net Weight	215kg
Digital Temperature Display	Yes
Digital Temperature Control	Yes
High/Low Temperature Alarm	Yes
Min/Max Temperature Recording	Yes
Door Open Alarm	Yes
Access Porthole	Optional Extra
Fan Assisted Cooling	Yes
Internal Light	Yes
Lockable Castors	Yes
Remote Alarm Contacts	Optional Extra
Auto Defrost	Yes
Free UKAS accredited temperature calibration Certificate	Included
Power Consumption Kwh/24hrs	3.6
Eco Friendly Refrigerant	Yes
Ambient Operating Temperature	16°C to 32°C

Features

- Digital temperature display
- Door open alarm
- Adjustable shelving
- Independent door locks
- Eco friendly refrigerant
- Double glazed insulating glass doors
- Free calibration certificate from our UKAS accredited calibration laboratory #8898
- High/low temperature alarm
- Automatic defrost
- Lockable castors
- LED energy saving Internal light
- Self closing doors
- Available as a solid door model see RPFR44043
- 2 year parts and labour warranty (UK only)

Models

	Model	Gross Capacity	Door Style	Lock Type
	Advanced Pharmacy Fridge RPF21043	650 litres	Solid	Key lock
	Advanced Pharmacy Fridge RPF21043	650 litres	Glass	Key lock
	Advanced Pharmacy Fridge RPF44043	1300 litres	Solid	2 key locks
	Advanced Pharmacy Fridge RPF44043	1300 litres	Glass	2 key locks



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.