



## RLVL03204

### Basic Laboratory Freezer



Constructed from enamelled steel, this freezer from the **Labcold** Basic Range is designed to provide basic storage at a lower price while still providing the confidence of the **Labcold** brand and support.

The freezer has two clear plastic drawers, one which is the interior width of the freezer and a half drawer at the bottom. There is a storage compartment that runs the whole depth of the freezer with a spring loaded flap to keep contents safe at the top. The adjustable feet ensures the freezer can be stabilised and manual defrost provides you with greater operator control over your temperature sensitive contents.

To put you in control, the temperature can be adjusted with an easy dial located on the right hand side of the interior, however, Labcold will by default set the set point to -20°C.

Please note, this freezer is not sparkfree and can not have a lock or digital temperature display fitted to it.

This Basic Freezer is ideal for economic cold storage of items that are not volatile and where precise storage temperature is not product critical. As you would expect from Labcold, this freezer comes backed with a full one years parts and labour warranty.\*



| Part Number                         | RLVL03204                                      |
|-------------------------------------|------------------------------------------------|
| Capacity (Litres)                   | 61                                             |
| Style                               | Underbench                                     |
| Exterior Dimensions (HxWxD)         | 849 x 475 x 448mm                              |
| Energy Consumption                  | 0.54 kWh/24h                                   |
| Noise level                         | 39dB(A)                                        |
| Temperature Control                 | ✓                                              |
| Operating Temperature               | -20°C                                          |
| Temperature Set Point Range         | -18°C to -23°C                                 |
| Sparkfree                           | ✗                                              |
| Interior Construction               | Moulded Plastic                                |
| Exterior Construction               | Enamelled Steel                                |
| Drawers                             | 2 + 1 flap                                     |
| Door Hinging                        | Right<br>(can be reversed at time of ordering) |
| Door Style                          | Solid                                          |
| Feet                                | 4 adjustable                                   |
| Manual Defrost                      | ✓                                              |
| GWP Friendly Refrigerant            | ✓                                              |
| Ammonia Free                        | ✓                                              |
| 1 year UK Parts and Labour Warranty | ✓                                              |

\* Labour UK only





# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.