

[Click Below](#)

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers

Rotators

MX-RD-Pro/MX-RD-E/MX-RL-Pro/MX-RL-E

DVAB

Your Accessible Solution

Very Flexible



MX-RD-Pro



MX-RD-E



MX-RL-Pro



MX-RL-E

Advantages for the user	Features
Different models	The available 4 types of rotators are round and linear, digital and economical
One machine double uses	The main body of MX-RD-Pro and MX-RL-Pro can be exchanged each other
Wide application	A wide variety of accessories are available

	MX-RD-E	MX-RD-Pro	MX-RL-E	MX-RL-Pro
Arguments				
Performance W/rpm	Two versions: for 10-70 (or 0-80 rpm s. below)			
Tilt Angle	Adjustable mixing angle between 0—90°		Mixing angle adjustment is not possible	
Comfort	Many accessories, up to 80 rpm	Speed control and optional timer, many accessories, LCD	Many accessories, up to 80 rpm	Speed control and optional timer, many accessories, LCD

Click Below

For General Use

Manual pipettes

Pip. Filler/ Stepper

Electronic Dispenser

Magnetic Stirrers

Shakers

For Analytical Labs

Electronic Titrator

Manual Dispenser

Rotary Evaporator

Hotplate

Overhead Stirrers

For Life Science Labs

Electronic Pipettes

SafeVac

Centrifuges-LS

Mixers / Vortex

Thermo Mixers

For Clinical Labs

Centrifuges-Clinical

Rollers and Rotators

Thermo Cyclers

Spectrophotometers

Rotators

MX-RD-Pro/MX-RD-E/MX-RL-Pro/MX-RL-E

DVAB

Your Accessible Solution

MX-RD-Pro

Technical data

Motor type	DC motor
Tilt angle	0-90°
Speed range	10-70rpm
Speed display	LCD
Timer display	LCD
Timer	1-1199min
Operation mode	Continuous / timer
Voltage	100-240V, 50/60Hz
Power	40W
Dimensions [D × W × H]	280 × 210 × 300mm
Permissible ambient temperature and humidity	5~40°C, 80%
Protection class	IP21



MX-RD-E

Technical data

Motor type	DC motor
Tilt angle	0-90°
Speed range	0-80rpm
Speed display	Scale
Timer display	Scale
Operation mode	Continuous
Voltage	100-240V, 50/60Hz
Power	20W
Dimensions [D×W×H]	300 × 220 × 310mm
Permissible ambient temperature and humidity	5~40°C, 80%
Protection class	IP21



Click Below

For General Use

Manual pipettes
Pip. Filler/ Stepper
Electronic Dispenser
Magnetic Stirrers
Shakers

For Analytical Labs

Electronic Titrator
Manual Dispenser
Rotary Evaporator
Hotplate
Overhead Stirrers

For Life Science Labs

Electronic Pipettes
SafeVac
Centrifuges-LS
Mixers / Vortex
Thermo Mixers

For Clinical Labs

Centrifuges-Clinical
Rollers and Rotators
Thermo Cyclers
Spectrophotometers

Rotators

MX-RD-Pro/MX-RD-E/MX-RL-Pro/MX-RL-E

DLAB

Your Accessible Solution

Varieties of accessories



Disk Accessories Type-1

Used for MX-RD-Pro
for 1.5mL x 60 centrifuge tubes holder

Cat. No.

18900160



Disk Accessories Type-2

Used for MX-RD-Pro
for 15mL x 16 centrifuge tubes holder

Cat. No.

18900161



Disk Accessories Type-3

Used for MX-RD-Pro
for 50mL x 8 centrifuge tubes holder

Cat. No.

18900162



Disk support rods

4 pcs in a group
use with disk accessories

Cat. No.

18900140



Rotisserie Accessory Type-1

Used for MX-RL-Pro, for 1.5mL x 48
centrifuge tubes held horizontally

Cat. No.

18900142



Rotisserie Accessory Type-2

Used for MX-RL-Pro, for 15mL x 24
centrifuge tubes held horizontally

Cat. No.

18900143



Rotisserie Accessory Type-3

Used for MX-RL-Pro, for 50mL x 24
centrifuge tubes held horizontally

Cat. No.

18900144



Rotisserie Accessory Type-4

Used for MX-RL-Pro, for 1.5mL x 32
centrifuge tubes held vertically

Cat. No.

18900145



Rotisserie Accessory Type-5

Used for MX-RL-Pro, for 15mL x 16
centrifuge tubes held vertically

Cat. No.

18900146



Rotisserie Accessory Type-6

Used for MX-RL-Pro, for 50mL x 16
centrifuge tubes held vertically

Cat. No.

18900147



Plastic Clips

Used with MX-RL-E/MX-RD-E
Disk/Rotisserie accessory, for
50mL, 15mL and 1.5mL
centrifuge tubes



Click Below

For General Use

Manual pipettes

Pip. Filler/ Stepper

Electronic Dispenser

Magnetic Stirrers

Shakers

For Analytical Labs

Electronic Titrator

Manual Dispenser

Rotary Evaporator

Hotplate

Overhead Stirrers

For Life Science Labs

Electronic Pipettes

SafeVac

Centrifuges-LS

Mixers / Vortex

Thermo Mixers

For Clinical Labs

Centrifuges-Clinical

Rollers and Rotators

Thermo Cyclers

Spectrophotometers

Rotators and Rollers

Ordering Information



Your Accessible Solution

Cat. No.	P/N	Descriptions
8031312000	MX-T6-S	Classic Tube Roller, variable speed, 6 rollers
8031312100	MX-T6-Pro	LCD Digital Tube Roller, variable speed, 6 rollers
8031412100	MX-RD-Pro	LCD Digital Rotator(including a 50ml disk accessory - 18900162), used with disk accessories,variable speed, adjustable mixing angle
8031402100	MX-RL-Pro	LCD Digital Rotator (including a 50ml rotisserie accessory - 18900147),used with rotisserie accessories,variable speed, adjustable mixing angle
8031412000	MX-RD-E	Rotator(including accessory-18900327/8/9), used with disk accessories,variable speed,
8031402000	MX-RL-E	Rotator(including accessory-18900327/8/9), used with rotisserie accessories,variable speed,
Accessories of Rotator		
Cat. No.	P/N	Descriptions
18900160		Disk accessory, for 1.5ml x 60 centrifuge tubes holder, use with MX-RD-Pro
18900161		Disk accessory, 15ml x 16 centrifuge tubes holder, use with MX-RD-Pro
18900162		Disk accessory, 50ml x 8 centrifuge tubes holder, use with MX-RD-Pro
18900140		Disk support rods, 4 pcs, use with disk accessories
18900142		Rotisserie accessory, for 1.5ml x 48 centrifuge tubes held horizontally, use with MX-RL-Pro
18900143		Rotisserie accessory, for 15ml x 24 centrifuge tubes held horizontally, use with MX-RL-Pro
18900144		Rotisserie accessory, for 50ml x 24 centrifuge tubes held horizontally, use with MX-RL-Pro
18900145		Rotisserie accessory, for 1.5ml x 32 centrifuge tubes held vertically, use with MX-RL-Pro
18900146		Rotisserie accessory, for 15ml x 16 centrifuge tubes held vertically, use with MX-RL-Pro
18900147		Rotisserie accessory, for 50ml x 16 centrifuge tubes held vertically, use with MX-RL-Pro
18200995		Disk and Rotisserie accessory, for 1.5mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E POM
18200996		Disk and Rotisserie accessory, for 15mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E POM
18200997		Disk and Rotisserie accessory, for 50mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E POM
18900327		Disk and Rotisserie accessory, for 1.5mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC
18900328		Disk and Rotisserie accessory, for 15mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC
18900329		Disk and Rotisserie accessory, for 50mlx 12 centrifuge tubes, use with MX-RL-E and MX-RD-E PC



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.