

# Chest freezer ULT C75



- Low cost of ownership
- Limited maintenance and on-site service
- Easily installable & movable

- Efficient cyclopentane insulation
- Eco-friendly, HCFC-free Nature R refrigerant
- ISO 9001, 14001 & 45001 Certified

## Specifications

### Technical

|                        |                  |
|------------------------|------------------|
| Gross volume           | 74 litres        |
| Net volume             | 71 litres        |
| Temp. range            | -60°C to -86°C   |
| Temp. probes           | 1/PT 1000        |
| Compressor             | Embraco          |
| Refrigerant            | Nature R         |
| Noise level            | 55 dB            |
| Power consumption 24/h | 5 kWh            |
| Power consumption      | 1.825 kWh / year |
| Voltage                | 230V             |
| Frequency              | 50 Hz            |

### Features

|                                             |        |
|---------------------------------------------|--------|
| Audible & visual alarm                      | YES    |
| Castors                                     | 4      |
| Integrated lock                             | YES    |
| Inner lids                                  | 1      |
| Defrost                                     | Manual |
| Battery Backup                              | 48h    |
| USB                                         | NO     |
| Data connection                             | NO     |
| Digital controller with display             | XR30CX |
| Data logging                                | NO     |
| High/Low temperature alarm                  | YES    |
| 2 portholes for external probe & CO2 backup | Ø20 mm |
| Prepared for external logging equipment     | YES    |
| Connection for remote alarm                 | YES    |

### Material

|       |                             |
|-------|-----------------------------|
| Inner | Precoated white steel       |
| Outer | Galvanized pre-coated steel |

### Dimensions

|                                  |        |
|----------------------------------|--------|
| Height (with castors)            | 870 mm |
| Width                            | 555 mm |
| Depth (incl. condenser)          | 580 mm |
| Depth (incl. condenser & handle) | 630 mm |
| Net weight                       | 45 kg  |

### Internal dimensions

|                   |        |
|-------------------|--------|
| Height            | 465 mm |
| Width             | 387 mm |
| Depth             | 387 mm |
| Insulation        | 80 mm  |
| Insulation type   | PUR    |
| Internal humidity | <5 %   |

## ULT C75



Product no.

N113001

## Storage capacity

| Box height | No. of boxes | No. of racks | Total no. of boxes | Product no. |
|------------|--------------|--------------|--------------------|-------------|
| 50 mm      | 8            | 5            | 40                 | N311210     |
| 75 mm      | 5            | 5            | 25                 | N311310     |
| 100 mm     | 4            | 5            | 20                 | N311410     |



# WolfLabs

**Pricing on any accessories shown can be found by keying the part number into the search box on our website.**

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

**[www.wolflabs.co.uk](http://www.wolflabs.co.uk)**

**Tel : 01759 301142**

**Fax : 01759 301143**

**[sales@wolflabs.co.uk](mailto:sales@wolflabs.co.uk)**

Please contact us if this literature doesn't answer all your questions.