

LAB

Seat or stool options for use in laboratory settings that are easily cleaned for maintaining hygiene



5 YEAR GUARANTEE FOR 8 HOUR USE



BACK ANGLE ADJUSTMENT

Lift the back right lever to free the tilt angle of chair back to your preferred position. To lock in position push the lever back down

SEAT HEIGHT ADJUSTMENT

Lift the front right lever to operate the gas lift. This will increase the chair height. Lifting whilst seated will lower chair height

DIMENSIONS

	PLAB(1) Black polyurethane lab chair	PLAB(2) Black polyurethane lab chair	PLAB(1)D Black polyurethane lab chair with draughtsman base	PLAB STOOL Black polyurethane lab stool	PLAB SMALL STOOL Black polyurethane small lab stool
SEAT HEIGHT	430-570mm	430-570mm	540-800mm	430-570mm	410-540mm
SEAT WIDTH	460mm	460mm	460mm	450mm	350mm
SEAT DEPTH	430mm	430mm	430mm	410mm	350mm
OVERALL HEIGHT	830-960mm	830-960mm	940-1200mm	530-670mm	410-540mm
OVERALL WIDTH	630mm	630mm	630mm	630mm	630mm
OVERALL DEPTH	630mm	630mm	630mm	630mm	630mm



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.