

1. Part Number: 5061A-250

2. Photo:



Closed Position

Note: Telescopic Rod is sold separately.



Open Position

Note: Telescopic Rod is sold separately.

3. Product Information:

Description	Snapper 250ml – 316Lss
Materials:	Half Spheres: 316L stainless steel Frame: 316L stainless steel Nuts & Bolts: 316L stainless steel Springs: 304 stainless steel
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	250ml
Max. Width when fully open:	275mm
Length:	195mm
Nominal Weight:	630g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

Note: Suitable for use with telescopic rods (Type 5061T)

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 5-Dec-24



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.