



Stomacher® 80 Biomaster series



Small tissue processing solutions for clinical and life science microbiology

The **Stomacher® 80 Biomaster** and **Stomacher® 80 microBiomaster** provide small tissue processing solutions for clinical, life science and microbiology applications. Capable of processing samples from 250ul to 80ml, the Biomaster range caters for the increased demand for smaller samples without reducing sample quality.

The **Stomacher® 80 Biomaster** (5ml to 80ml sample capacity), and **Stomacher® Bags** allow samples to be processed safely and without risk of cross contamination while maximising visible cell recovery. Tens of thousands of referenced methods are available and now with the enhanced

sample capacity range of the microBiomaster (250ul to 80ml sample capacity) the **Stomacher® 80** is the World's most efficient and versatile sample preparation system for life science applications.

Choose the **Stomacher® 80 microBiomaster** (250ul to 80ml sample capacity) for increased sampling volume range and adjustable process times & paddle speeds.

Why the Stomacher® 80?

Adjustable process times
& paddle speeds

Proven enhancement in
cell recovery & viability

No risk of cross
contamination

Process samples
from 250ul - 80ml

Cost effective sample
preparation



Stomacher® 80
Biomaster



Stomacher® 80
microBiomaster


Stomacher® 80 Biomaster



The **Stomacher® Bag** Range has been designed to fulfil all applications from blending, straining and storage including bags designed specifically for life science applications.

Enhanced **Stomacher® 80 Bag** range features closure and strainer bag formats and double bags for added biosafety.

Seward Bag racks provide the perfect storage for **Stomacher® Bag** sample handling, helping to provide a clean and productive working environment.



Stomacher® 80

Product Code	Description	Voltage	Pack Dimensions (LxWxH)	Weight (kg)
0080/000/EU	Stomacher® 80 Biomaster	240V	49.5x33.5x23.5cm	18.0
0080/000/AJ	Stomacher® 80 Biomaster	110V	49.5x33.5x23.5cm	18.5
0080/000/EU/M	Stomacher® 80 microBiomaster	240	49.5x33.5x23.5cm	18.0
0080/000/AJ/M	Stomacher® 80 microBiomaster	110V	49.5x33.5x23.5cm	18.5

80 Bags

Product Code	Bag Type	Bag Size	Sample Size	Pack Quantity
BA6040	Standard Bag	105 x 155mm	250ul-80ml	2000 (sachets of 50)
BA6040CLR	Closure Bag	105 x 155mm	5-80ml	200 (sachets of 10)
BA6040/STR	Strainer Bag	105 x 155mm	5-80ml	200 (sachets of 10)
BA6040/STR/DBL	Strainer Double Bag	135 x 180mm	5-80ml	200 (sachets of 10)
BA6040/CLR/STR	Strainer Closure Bag	105 x 155mm	5-80ml	200 (sachets of 10)

Accessories

Product Code	Description	Capacity	Dimensions	Weight
BA6098	Stomacher® 80 Bag Opener	1 Place	7 x 7 x 14cm	0.1kg
BA6090	Stomacher® 80 Bag Rack	10 Place	30 x 11 x 12cm	0.9kg
BA6051	Simplettes® Disposable straws	200 (8 sachets of 25 each)	24 x 18 x 18cm	3.7kg



WolfLabs

Pricing on any accessories shown can be found by keying the part number into the search box on our website.

The specifications listed in this brochure are subject to change by the manufacturer and therefore cannot be guaranteed to be correct. If there are aspects of the specification that must be guaranteed, please provide these to our sales team so that details can be confirmed.

www.wolflabs.co.uk

Tel : 01759 301142

Fax : 01759 301143

sales@wolflabs.co.uk

Please contact us if this literature doesn't answer all your questions.